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INSTRUCTIONS FOR USE DIAMIR dental diamond burs CE 0051

Field of use

- Preparation of the cavity (preservative)
- Preparation of the stumps (prosthetic)
- Finishing of the cavity edges and repair work
- Polishing of the repairs

Instructions for use

To get the best from **DIAMIR** instruments, the instructions listed herewith should be followed carefully:

- Use suitable cooling methods to guarantee correct cooling of the parts subject to treatment.
- Use aspiration techniques which do not compromise cooling of the parts involved in the treatment.
- Use only hand pieces in perfect condition, as well as cooling nozzles in perfect operation.
- Use only tools which are in perfect condition, not only in terms of the diamond edge but also as regards concentricity.
- Insert the instruments very carefully and without forcing.
- Before using, check that the tool shank locks perfectly with the hand piece.
- Bring the tool to the necessary speed before approaching the part to be treated.
- Use a pressure force between 0.3 - 2 N (approximately 30 - 200 g) Under no circumstance must greater pressures be applied.
- Avoid using the tools as levers and/or at incorrect angulations.
- Respect the speeds shown in the attached table.

Advice and Recommendations on Hygiene

Before using on the patient, make sure that the tools are disinfected, cleaned and sterilised. Obviously this must be repeated before each use.

The tools must be disinfected, cleaned and sterilised according to proven medical protocols.

It is recommended to sterilized the tools by steam, Autoclave type B, conform to EN 13060, with minimum sterilization cycle of 5 minutes at 134° C degrees.

In any case, use the sterilization devices (materials and / or equipment) respecting the instructions provided by the manufacturer. The sterility of the instruments is the user's responsibility.

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